

volume to resin volume is approximately eight parts load volume to one part C-18 resin;  
and

formulating the collected protein solution for final use.

a<sup>3</sup> 3. (Amended) The method of claim 1, wherein the ethanol solution is  
comprised of about 19% ethanol and about 81% purified water adjusted to a pH of 5.7 to  
5.8.

IN THE ABSTRACT

Please enter the following substitute paragraph for the abstract at page 15,  
line 1 through line 7:

a<sup>4</sup> The method for immune serum globulin purification relates to the  
purification of immune globulins from blood plasma with a high degree of efficiency and  
a high rate of recovery. The immune globulin source is Cohn's fraction I+II+III or II+III  
prepared from plasma or plasma intermediates by precipitation of the paste at pH 5.7 to  
5.8 in the presence of 20% ethanol and 80% purified water. A glycine extraction is  
followed by an anion exchange chromatography column step to achieve a significantly  
high yield and high purity of the concentrated protein.